

Application by Health Care Provider for Enrolment in the Vaccination Subsidy Scheme (VSS)

Record of the Vaccine Storage Refrigerator


Doctors enrolled in VSS have to use Purpose-built vaccine refrigerators (PBVR) for vaccine storage by 2023/24 vaccination season. For further details regarding vaccine storage under VSS, please refer to Chapter 6 of the VSS Doctor's Guide – Requirements on vaccine storage and handling under VSS (https://www.chp.gov.hk/files/pdf/vssdg_ch6_vaccine_storage_and_handling.pdf).

For record and follow-up purpose, please complete this form, supplemented with:-

- (1) One photo of the interior of the refrigerator **currently** used for vaccine storage at your clinic; and
- (2) One photo of the temperature monitoring device AND temperature record

Should you have any enquiries, please contact Ms. TAM at 3975 4848 or Ms. CHAN at 3975 4846.

(1) Photo of the interior of the refrigerator and Max-min thermometer



Fridge Brand Name: _____

Fridge Model No.: _____

Doctor's Name: _____

Practice Name and No.: _____

SPID: _____

Contact Phone number: _____

Signature: _____

Date: _____

(2) Photo or copy of the temperature record kept in the clinic

貯存疫苗的雪櫃溫度檢查表

An Example only

1. 所有疫苗(季節性流感疫苗及肺炎球菌疫苗)，須保存於攝氏 +2℃ 至 +8℃ 雪櫃內備用。
2. 須每日檢查 三次（上午、中午及下午 各一次） 雪櫃內的 現時溫度、雪櫃最高及最低溫度，並記錄在本表格上。請於每次記錄雪櫃溫度後，重置最高最底溫度計。
3. 請保留此記錄至少一年，以便有需要時作參考。

**註: 如雪櫃溫度 低於攝氏+2 度 或 高於攝氏+8 度

1. 請暫勿使用受影響的疫苗，並應將疫苗立即存放於攝氏+2 度至+8 度的雪櫃。
2. 請通知診所醫生及向疫苗供應商查詢
3. 在' 溫度超溫事故報告表' 填寫有關事故並傳真至項目管理及疫苗計劃科

日期	檢查 雪櫃 時間	雪櫃內溫度			記錄人員			備註
		現時	最高	最低	姓名	職位	簽署	
15-08-2022	09:00	5.3	6	4.1				
15-08-2022	13:00	5.4	6	4.1				
15-08-2022	17:00	5.4	6	4.1				